



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® screw caps ND 13 with hole, septum
natural rubber red-or./TEF transp., 100
pcs./pack**

**€19.80
plus VAT &
Shipping**

Product Images



Description

Black screw caps ND13 made of PP, with thread 13-425.

With hole 8.5 mm.

With septum natural rubber red-orange/TEF transparent

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness 1.3 mm

Additional Information

No.	7-0770
Manufacturer (Brand)	neochrom
EAN	4058072030193
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	Reference to non-verified autosampler compatibility in shopping cart suitable for autosampler suitable for multiple injections cap with septum cap autoclavable
DM Hole	8.5 mm
Thread	13-425 (US)
Quantity components	two
Hardness	60 shore A
Material gasket / septum	Natural rubber
Texture	corrugated
Closure / lid properties	with hole with mounted septum
DM Closure	13 mm
Thickness septum	1.3 mm
Color closure/ lid	black
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

