



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neochrom® screw caps ND 13 with hole, septum silicone cream/PTFE red, 100 pcs./pack

**€36.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Black screw caps ND13 made of PP, with thread 13-425.  
With hole 8.5 mm.

With septum silicone cream/PTFE red  
Temperature resistant from -60°C to +200°C  
For single use injections  
Thickness 1.5 mm

## Additional Information

No.	7-0772
Manufacturer (Brand)	neochrom
EAN	4058072030216
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	Note on non-verified autosampler compatibility in shopping cart suitable for autosampler cap with septum cap autoclavable
Quantity components	two
Texture	corrugated
DM Hole	8.5 mm
DM Closure	13 mm
Thickness septum	1.5 mm
Color closure/ lid	black
Thread	13-425 (US)
Hardness	55 shore A
Material gasket / septum	Silicone Polytetrafluoroethylene (PTFE)
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
Closure / lid properties	with hole with mounted septum
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

