



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neochrom® UltraBond screw caps ND24 w. Hole, silicone natural/PTFE beige, 100 pcs./pack

**€64.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Cap and septum form an inseparable unit in the <sup>UltraBond™</sup> closure caps, so that the septum cannot be poked into the bottle with a blunt needle.

Screw caps made of white PP with thread 24-400, with hole 15 mm Ø or closed.

The septum silicone/PTFE in EPA quality is temperature resistant from -60°C to +200°C and suitable for single injections.

Thickness 3.2 mm

Hardness 45° shore A

Order No. 7-0864

Screw cap with hole Ø 15mm

## Additional Information

No.	7-0864
Manufacturer (Brand)	neochrom
EAN	4058072037321
Transport temperature	Room temperature
Color	White
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	Reference to non-verified autosampler compatibility in shopping cart suitable for autosampler cap with septum cap autoclavable septum and cap one unit EPA certified
DM Hole	12.5 mm
Thread	24-400 (US)
Quantity components	two
Hardness	45 shore A
Material gasket / septum	Silicone Polytetrafluoroethylene (PTFE)
Texture	corrugated
Closure / lid properties	with hole with mounted septum
DM Closure	24 mm
Thickness septum	3.2 mm
Color closure/ lid	white
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

