



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® UltraBond screw caps ND24 closed, silicone natural/PTFE beige, 100 pcs./pack

€69.50
plus VAT &
Shipping

Product Images



Description

Cap and septum form an inseparable unit in the ^{UltraBond™} closure caps, so that the septum cannot be poked into the bottle with a blunt needle.

Screw caps made of white PP with thread 24-400, with hole 15 mm Ø or closed.

The septum silicone/PTFE in EPA quality is temperature resistant from -60°C to +200°C and suitable for single injections.

Thickness 3.2 mm

Hardness 45° shore A

Order No. 7-0865

Screw cap closed

Additional Information

No.	7-0865
Manufacturer (Brand)	neochrom
EAN	4058072037338
Transport temperature	Room temperature
Color	White
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	Reference to non-tested autosampler compatibility in shopping cart Cap with septum Cap autoclavable Septum and cap one unit EPA certified
Thread	24-400 (US)
Quantity components	two
Hardness	45 shore A
Material gasket / septum	Silicone Polytetrafluoroethylene (PTFE)
Texture	corrugated
Closure / lid properties	closed with mounted septum
DM Closure	24 mm
Thickness septum	3.2 mm
Color closure/ lid	white
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

