



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neochrom® screw caps ND8 closed, butyl red/PTFE grey, 100 pcs./pack

**€16.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Black screw caps ND8 made of PP, with thread 8-425, closed.

With septum Butyl red/PTFE grey

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness 1.3 mm

## Additional Information

No.	7-0872
Manufacturer (Brand)	neochrom
EAN	4058072037383
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	Closure with septum Closure autoclavable
Thread	08-425 (US)
Quantity components	two
Hardness	55 shore A
Material gasket / septum	Butyl rubber (IIR) polytetrafluoroethylene (PTFE)
Texture	corrugated
Closure / lid properties	closed with mounted septum
Color closure/ lid	black
DM Closure	8 mm
Thickness septum	1.3 mm
Spare part	2-7203
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

