



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® screw caps ND8 closed, rubber red- or./TEF transp., 100 pcs./pack

€14.90
plus VAT &
Shipping

Product Images



Description

Black screw caps ND8 made of PP, with thread 8-425, closed.

With septum Natural rubber red-orange/TEF transparent

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness 1.3 mm

Additional Information

No.	7-0873
Manufacturer (Brand)	neochrom
EAN	4058072037390
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	Closure with septum Closure autoclavable
Quantity components	two
Texture	corrugated
DM Closure	8 mm
Thickness septum	1.3 mm
Color closure/ lid	black
Thread	08-425 (US)
Hardness	60 shore A
Material gasket / septum	Natural rubber Polytetrafluoroethylene (PTFE)
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
Closure / lid properties	closed with mounted septum
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

