



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® Screw Caps ND8 w. Hole, PP,
Septum Red Rubber/PTFE beige, 100 pcs./pack**

**€16.50
plus VAT &
Shipping**

Product Images



Description

Black screw caps ND8 made of PP, with thread 8-425 and hole 5.5 mm.

1. with septum natural rubber/TEF and butyl/PTFE

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness 1.3 mm

2. with septum silicone/PTFE

Temperature resistant from -60°C to +200°C

For single injections

Thickness 1.0 mm

Additional Information

No.	EC-1300
Manufacturer (Brand)	neochrom
EAN	4058072085926
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	suitable for autosampler cap with septum cap autoclavable
DM Hole	5.5 mm
Quantity components	two
Hardness	45 shore A
Material gasket / septum	Polytetrafluoroethylene (PTFE) rubber
Texture	corrugated
Closure / lid properties	with hole with mounted septum
Color closure/ lid	black
Thickness septum	1 mm
TBST MAX Septum	200 °C
TBST MIN Septum	-40 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

