



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® Screw Caps ND8 w. Hole, PP, silicone
cream/PTFE red Ultra Clean, 100 pcs./pack**

**€30.50
plus VAT &
Shipping**

Product Images



Description

Black screw caps ND8 made of PP, with thread 8-425 and hole 5.5 mm.

1. with septum natural rubber/TEF and butyl/PTFE

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness 1.3 mm

2. with septum silicone/PTFE

Temperature resistant from -60°C to +200°C

For single injections

Thickness 1.0 mm

Additional Information

No.	EC-1301
Manufacturer (Brand)	neochrom
EAN	4058072085933
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	with high purity septum for critical analyses suitable for autosampler cap with septum cap autoclavable
Quantity components	three
Texture	corrugated
DM Hole	5.5 mm
Thickness septum	1.5 mm
Color closure/ lid	black
Hardness	45 shore A
Material gasket / septum	Silicone Polytetrafluoroethylene (PTFE)
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
Closure / lid properties	with hole with mounted septum
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

