



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® Screw cap ND9 w. Hole, blue,
septum silicone white/PTFE red, slotted (y) 100
pcs./pack**

**€45.90
plus VAT &
Shipping**

Product Images



Description

Colored screw caps ND9 made of PP, with short thread 9-425 and hole 6 mm.

1. with septum natural rubber/TEF and red rubber/PTFE

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness 1.0 mm

2. with septum silicone/PTFE

Temperature resistant from -60°C to +200°C

For single injections

Thickness 1.0 mm

3. with septum PTFE virginal

high temperature resistant

For single use injections

Thickness 0.2 mm

Additional Information

No.	EC-1305
Manufacturer (Brand)	neochrom
EAN	4058072085964
Transport temperature	Room temperature
Color	blue
Material	Polypropylene (PP)
Closures and lids features	with slotted septum for blunt. thick needles suitable for autosampler cap with septum cap autoclavable septum with low particle generation
Quantity components	three
Texture	corrugated
DM Hole	6 mm
DM Closure	9 mm
Thickness septum	1 mm
Color closure/ lid	blue
Thread	Short thread
Hardness	55 shore A
Material gasket / septum	Polytetrafluoroethylene (PTFE) Silicone
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
Closure / lid properties	with hole with mounted septum
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

