



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® Short Thread MS Cap, transparent, with diaphragm

€34.30
plus VAT &
Shipping

Product Images



Description

MS-Cap with short thread 9-425, non-bleeding one-component cap for LC/MS and GC/MS.
Transparent, with diaphragm.

Additional Information

No.	7-0714
Manufacturer (Brand)	neochrom
EAN	4058072029739
Transport temperature	Room temperature
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	Reference to untested autosampler compatibility in shopping cart suitable for mass spectrometry suitable for autosampler
DM Hole	5.5 mm
Thread	Short thread
Texture	corrugated
Closure / lid properties	pierceable
DM Closure	9 mm
Color closure/ lid	transparent (milky)
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

