



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz octagonal cover plug NS 12/21, solid design

**€94.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Octagonal glass lid stopper

According to DIN 12252 Form A and E, in solid design.

## Additional Information

---

No.	LZ-0257
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2256
EAN	4058072299996
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (crystal clear)
Material	Borosilicate glass
sterile	No
Type closure/ lid	NS plug
suitable for	Bottles with standard ground joint
Closures and lids features	with standard ground joint
Number of edges	8
Standard ground core	12/21
TBST MAX	200 °C
TBST MIN	-70 °C
standardized according to DIN standard (DIN)	DIN 12252

