



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz octagonal lid plug NS 45/40, semi-hollow version

€27.50
plus VAT &
Shipping

Product Images



Description

Octagonal glass lid stopper

According to DIN 12252 Form A and E, in semi-hollow design.

Additional Information

No.	LZ-0263
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2262
EAN	4058072300050
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (crystal clear)
Material	Borosilicate glass
sterile	No
Type closure/ lid	NS plug
suitable for	Bottles with standard ground joint
TBST MAX	200 °C
TBST MIN	-70 °C
Closures and lids features	with standard ground joint
Number of edges	8
Standard ground core	45/40
standardized according to DIN standard (DIN)	DIN 12252

