



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Starke hexagonal hollow plugs brown, NS 12/21, 5 pcs./pack

**Price on  
request  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

### **Conical hexagonal hollow glass stoppers for all vessels with normal ground joint.**

Fully resistant to temperature changes, very light.

Can also be used for light-sensitive media.

## Additional Information

---

No.	SE-0952
Manufacturer (Brand)	Starke
old neoLab article no.	2-4267
EAN	4066292206977
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	brown
Material	Borosilicate glass
sterile	No
Type closure/ lid	NS plug
suitable for	Bottles with standard ground joint
TBST MAX	200 °C
TBST MIN	-70 °C
Closures and lids features	with standard ground joint
Number of edges	6
Standard ground core	12/21

