



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Drosophila stoppers, suitable for 177
ml vials, white, foam, 100 pcs/pack**

**€9.90
plus VAT &
Shipping**

Product Images



Description

Plugs for Drosophila rearing.

Suitable plugs for the vessels or vials made of high-density foam. This allows foreign particles to be filtered and oxygen to diffuse in.

Additional Information

No.	4-0073
Manufacturer (Brand)	neoLab
EAN	4058072084974
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	38,5 mm
sterile	No
Form	cylindrical
Height	45 mm
TBST MAX	80 °C
TBST MIN	-30 °C

