




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz octagonal lid plug NS 24/29, PE

€11.30
plus VAT &
Shipping

Product Images



Description

Plastic **octagonal lid plug**

According to DIN 12254 and in addition to the DIN standard.

Material:PE; Color: white-red

Additional Information

No.	LZ-0269
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2268
EAN	4058072300111
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polyethylene (PE)
sterile	No
Type closure/ lid	NS plug
suitable for	Bottles with standard ground joint
TBST MAX	200 °C
TBST MIN	-70 °C
Closures and lids features	with standard ground joint
Number of edges	8
Standard ground core	24/29
standardized according to DIN standard (DIN)	DIN 12254

