



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® silicone sponge stopper white 12 x 8 x 20 mm, 10 pcs./pack

**€23.20**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The plugs have the chemical and thermal properties of silicone and are lightweight as a feather.

Suitable for all round necks.

The stoppers are as easy to drill as a cork stopper.

Temperature resistant from -60°C to +200°C, autoclavable, odorless and non-toxic.

## Additional Information

---

No.	2-4250
Manufacturer (Brand)	neoLab
EAN	4058072002787
Transport temperature	Room temperature
autoclavable at 121°C	Yes
sterile	No
Form	conical
Height	20 mm
Type closure/ lid	Plug
DM top	12 mm
DM below	8 mm
TBST MAX	200 °C
TBST MIN	-60 °C
for NS	44531

