



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® silicone sponge stopper yellow 25 x 13  
x 17.5 mm, 10 pcs./pack**

**€22.70  
plus VAT &  
Shipping**

## Product Images



## Description

---

The plugs have the chemical and thermal properties of silicone and are lightweight as a feather.

Suitable for all round necks.

The stoppers are as easy to drill as a cork stopper.

Temperature resistant from -60°C to +200°C, autoclavable, odorless and non-toxic.

## Additional Information

---

No.	2-4252
Manufacturer (Brand)	neoLab
EAN	4058072002800
Transport temperature	Room temperature
autoclavable at 121°C	Yes
sterile	No
Form	conical
Height	25 mm
Type closure/ lid	Plug
TBST MAX	200 °C
TBST MIN	-60 °C
DM top	17.5 mm
DM below	13 mm

