



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® silicone sponge stopper pink 35 x 28 x 35 mm, 5 pcs./pack

€39.80
plus VAT &
Shipping

Product Images



Description

The plugs have the chemical and thermal properties of silicone and are lightweight as a feather.

Suitable for all round necks.

The stoppers are as easy to drill as a cork stopper.

Temperature resistant from -60°C to +200°C, autoclavable, odorless and non-toxic.

Additional Information

No.	2-4256
Manufacturer (Brand)	neoLab
EAN	4058072002848
Transport temperature	Room temperature
autoclavable at 121°C	Yes
sterile	No
Form	conical
Height	35 mm
Type closure/ lid	Plug
TBST MAX	200 °C
TBST MIN	-60 °C
DM top	35 mm
DM below	28 mm

