



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® silicone sponge stopper yellow 40 x 33 x 40 mm, 5 pcs./pack

€50.30
plus VAT &
Shipping

Product Images



Description

The plugs have the chemical and thermal properties of silicone and are lightweight as a feather.

Suitable for all round necks.

The stoppers are as easy to drill as a cork stopper.

Temperature resistant from -60°C to +200°C, autoclavable, odorless and non-toxic.

Additional Information

| | |
|-----------------------|------------------|
| No. | 2-4257 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072002855 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| sterile | No |
| Form | conical |
| Height | 40 mm |
| Type closure/ lid | Plug |
| TBST MAX | 200 °C |
| TBST MIN | -60 °C |
| DM top | 40 mm |
| DM below | 33 mm |

