



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® silicone sponge stopper blue 48 x 44 x 52 mm, 1 piece

**€25.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The plugs have the chemical and thermal properties of silicone and are lightweight as a feather.

Suitable for all round necks.

The stoppers are as easy to drill as a cork stopper.

Temperature resistant from -60°C to +200°C, autoclavable, odorless and non-toxic.

## Additional Information

---

No.	2-4259
Manufacturer (Brand)	neoLab
EAN	4058072002879
Transport temperature	Room temperature
autoclavable at 121°C	Yes
sterile	No
Form	conical
Height	48 mm
Type closure/ lid	Plug
TBST MAX	200 °C
TBST MIN	-60 °C
DM top	52 mm
DM below	44 mm

