



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Steri Stopper for single use, for neck 28 mm

€65.70
plus VAT &
Shipping

Product Images



Description

Disposable sterile stopper made of cellulose for sterile closure of test tubes, Erlenmeyer flasks and culture flasks.

Autoclavable up to 200°C.

Air-permeable, bacteria-filtering.

Additional Information

| | |
|-----------------------|---------------|
| No. | 1-1357 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072032630 |
| autoclavable at 121°C | Yes |
| sterile | No |
| Form | conical |
| Height | 30 mm |
| DM top | 28 mm |
| DM below | 23 mm |
| TBST MAX | 200 °C |

