



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Deutsch & Neumann rubber stopper red 35 x 29 mm 30 mm high; PU: 5 pieces

**€13.30** As low as  
**€10.03**  
plus VAT &  
Shipping

### Product Images



## Description

---

Made of natural rubber Para red, according to DIN 12871. Low sulfur, floating, non-toxic. 40 Shore A. Packing size according to requirements.

## Additional Information

---

No.	NM-0037
Manufacturer (Brand)	Deutsch & Neumann
old neoLab article no.	1-1021
EAN	4058072336974
Transport temperature	Room temperature
autoclavable at 121°C	Yes
sterile	No
Form	conical
Height	30 mm
Type closure/ lid	Plug
DM top	35 mm
DM below	29 mm

