



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Deutsch & Neumann rubber stopper red 44 x 36 mm 40 mm high; PU: 5 pieces

€22.66 As low
as
€17.09
plus VAT &
Shipping

Product Images



Description

Made of natural rubber Para red, according to DIN 12871. Low sulfur, floating, non-toxic. 40 Shore A. Packing size according to requirements.

Additional Information

No.	NM-0039
Manufacturer (Brand)	Deutsch & Neumann
old neoLab article no.	1-1023
EAN	4058072336998
Transport temperature	Room temperature
autoclavable at 121°C	Yes
sterile	No
Form	conical
Height	40 mm
Type closure/ lid	Plug
DM top	44 mm
DM below	36 mm

