



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Deutsch & Neumann rubber stopper red 59.5 x 50.5 mm 45 mm high; PU: 5 pieces

**€38.51**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Made of natural rubber Para red, according to DIN 12871. Low sulfur, floating, non-toxic. 40 Shore A. Packing size according to requirements.

## Additional Information

---

No.	NM-0042
Manufacturer (Brand)	Deutsch & Neumann
old neoLab article no.	1-1026
EAN	4058072337025
Transport temperature	Room temperature
autoclavable at 121°C	Yes
sterile	No
Form	conical
Height	45 mm
Type closure/ lid	Plug
DM top	59.5 mm
DM below	50.5 mm

