




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Inverting stopper 10 mm \varnothing NS 10 made of silicone PU 10 pieces

€11.20
plus VAT &
Shipping

Product Images



Description

Silicone stopper for absolute tightness even with defective glass rim. Material removal via piercing. Tight fit by turning the rim upside down.

Additional Information

No.	1-1034
Manufacturer (Brand)	neoLab
EAN	4058072031510
autoclavable at 121°C	Yes
sterile	No
Form	conical
Type closure/ lid	Plug NS Plug
DM below	10 mm
Closure / lid properties	with invertible rim pierceable
for NS	43739

