



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Inverting stopper 13 mm dia. NS 12.5
made of silicone PU 10 pcs. No. 1-1036**

**€17.40
plus VAT &
Shipping**

Product Images



Description

Silicone stopper for absolute tightness even with defective glass rim. Material removal via piercing. Tight fit by turning the rim upside down.

Additional Information

No.	1-1035
Manufacturer (Brand)	neoLab
EAN	4058072031527
autoclavable at 121°C	Yes
sterile	No
Form	conical
Type closure/ lid	Plug NS Plug
DM below	13 mm
Closure / lid properties	with invertible rim pierceable
for NS	44531

