




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Inverting plug 15 mm ø NS 14.5 made  
of silicone PU 10 pcs.**

**€22.20  
plus VAT &  
Shipping**

## Product Images



## Description

---

Silicone stopper for absolute tightness even with defective glass rim. Material removal via piercing. Tight fit by turning the rim upside down.

## Additional Information

---

No.	1-1036
Manufacturer (Brand)	neoLab
EAN	4058072031534
autoclavable at 121°C	Yes
sterile	No
Form	conical
Type closure/ lid	Plug NS Plug
DM below	15 mm
Closure / lid properties	with invertible rim pierceable
for NS	14/23

