



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Inverting stopper 25 mm ø NS 24 made of silicone PU 10 pcs. No. 1-1038

**€55.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Silicone stopper for absolute tightness even with defective glass rim. Material removal via piercing. Tight fit by turning the rim upside down.

## Additional Information

---

No.	1-1037
Manufacturer (Brand)	neoLab
EAN	4058072031541
autoclavable at 121°C	Yes
sterile	No
Form	conical
Type closure/ lid	Plug NS Plug
DM below	25 mm
Closure / lid properties	with invertible rim pierceable
for NS	24/29

