



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® NS stopper made of PE 8-Kt plate f. NS 12/21 PU 10 pcs. E-1584

Special Price

**€2.10** was

**€8.50**

**plus VAT &  
Shipping**

### Product Images



## Description

---

Cone smooth,  
Octagonal lid plug

## Additional Information

---

|                            |                                    |
|----------------------------|------------------------------------|
| No.                        | E-1584                             |
| Manufacturer (Brand)       | neoLab                             |
| EAN                        | 4058072065768                      |
| Transport temperature      | Room temperature                   |
| autoclavable at 121°C      | Yes                                |
| Color                      | White                              |
| Material                   | Polyethylene (PE)                  |
| sterile                    | No                                 |
| Type closure/ lid          | NS plug                            |
| suitable for               | Bottles with standard ground joint |
| TBST MAX                   | 200 °C                             |
| TBST MIN                   | -70 °C                             |
| Closures and lids features | with standard ground joint         |
| Number of edges            | 8                                  |
| Standard ground core       | 12/21                              |

