



neoLab Migge GmbH

Rischerstr. 7-9

69123 Heidelberg

Deutschland

+49 (0)6221 / 8442-44

<https://www.neolab.de>

Umsatzsteuer-Identifikationsnummer:

DE 143 450 657



neoLab® NS stopper made of PE 8-Kt plate for NS 14/23, PU 10 pieces

€8.50

plus VAT & Shipping

Product Images



Description

Cone smooth,
Octagonal lid plug

Additional Information

No.	E-1585
Manufacturer (Brand)	neoLab
EAN	4058072065775
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polyethylene (PE)
sterile	No
Type closure/ lid	NS plug
suitable for	Bottles with standard ground joint
TBST MAX	200 °C
TBST MIN	-70 °C
Closures and lids features	with standard ground joint
Number of edges	8
Standard ground core	14/23

