



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® NS stopper made of PE 8-Kt plate for NS 29/32 PU 10 pieces

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Cone smooth,
Octagonal lid plug

Additional Information

No.	E-1588
Manufacturer (Brand)	neoLab
EAN	4058072065805
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polyethylene (PE)
sterile	No
Type closure/ lid	NS plug
suitable for	Bottles with standard ground joint
TBST MAX	200 °C
TBST MIN	-70 °C
Closures and lids features	with standard ground joint
Number of edges	8
Standard ground core	29/32

