



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand Lid PE, f. tubes 1.2 ml, PP, strips of 8, non-sterile

**€37.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The autoclavable tubes can be used in many ways in the laboratory, whether for PCR, EIA, storage or transport.

The vials are freezable down to  $-80^{\circ}\text{C}$ .

The lids made of PE are not suitable for autoclaving.

Dimensions of the tubes:

- Outer diameter: 8.8 mm

- Height: 45 mm

## Additional Information

---

No.	BN-2286
Manufacturer (Brand)	Brand
old neoLab article no.	7-5243
EAN	4058072214951
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
sterile	No
Analytical method	Standard PCR
Number of tubes	8
Color closure/ lid	transparent (crystal clear)
Material closure/ lid	Polyethylene (PE)
Closure / lid properties	closed

