



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand lid with condensation rings for plates,
96-well with transparent bottom, VE=100, 5
pcs./btl.**

€78.45 As low
as
€74.40
plus VAT &
Shipping

Product Images



Description

To protect samples, reaction mixtures or cell cultures from contamination and evaporation, matching lids are available for all ^{BRANDplates®}. The lids are made of crystal-clear polystyrene. They are optimally matched to the corresponding ^{BRANDplates®} so that plates with lids can also be stacked easily. The cut-off corners provide the orientation of the lid.

Additional Information

No.	BN-2382
Manufacturer (Brand)	Brand
old neoLab article no.	7-5388
EAN	4058072215897
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Material	Polystyrene (PS)
sterile	No
Wide	80 mm
Length	80 mm
suitable for	Microtiter plates
TBST MAX	100 °C
Dimensions	80 x 130 mm

