



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand lid without condensation rings for plates,
96-well, standard, VE=100, 5 pcs./bag.**

€71.25 As low as
€67.60
**plus VAT &
Shipping**

Product Images



Description

To protect samples, reaction mixtures or cell cultures from contamination and evaporation, matching lids are available for all BRANDplates®. The lids are made of crystal-clear polystyrene. They are optimally matched to the corresponding BRANDplates® so that plates with lids can also be stacked easily. The cut-off corners provide the orientation of the lid.

Additional Information

No.	BN-2379
Manufacturer (Brand)	Brand
old neoLab article no.	7-5385
EAN	4058072215866
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Material	Polystyrene (PS)
sterile	No
Wide	80 mm
Length	80 mm
suitable for	Microtiter plates
TBST MAX	100 °C
Dimensions	80 x 130 mm

