



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Fire lid without condensation rings for plates 96-well with transparent bottom and all 384- well plates

**€69.35** As low  
as  
**€65.95**  
plus VAT &  
Shipping

### Product Images



## Description

---

To protect samples, reaction mixtures or cell cultures from contamination and evaporation, matching lids are available for all BRANDplates®. The lids are made of crystal-clear polystyrene. They are optimally matched to the corresponding BRANDplates® so that plates with lids can also be stacked easily. The cut-off corners provide the orientation of the lid.

## Additional Information

---

No.	BN-2380
Manufacturer (Brand)	Brand
old neoLab article no.	7-5386
EAN	4058072215873
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Material	Polystyrene (PS)
sterile	No
Wide	80 mm
Length	80 mm
suitable for	Microtiter plates
TBST MAX	100 °C
Dimensions	80 x 130 mm

