



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Fire sealing foil for plates, PP, DMSO-resistant,
self-adhesive, enables visual control**

**€98.70
plus VAT &
Shipping**

Product Images



Description

Self-adhesive disposable films can be used to securely cover and seal samples.

PP sealing film,

, which covers a temperature range of -40 to +120°C, is particularly suitable for storage for PCR and ELISA.

Aluminum sealing films are suitable for cold storage due to their good temperature resistance from -80 to +120°C.

Additional Information

No.	BN-1529
Manufacturer (Brand)	Brand
old neolab article no.	7-5235
EAN	4058072208868
Transport temperature	Room temperature
Color	White
Material	Polypropylene (PP)
Temp MAX	120 °C
Temp MIN	-40 °C
sterile	No
Format/ Packing	bundled packed

