



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Fire sealing foil for plates, aluminum, self-adhesive, for cold storage


Special Price

**€99.70** was

**€136.60**

**plus VAT &  
Shipping**

### Product Images



## Description

---

Self-adhesive disposable films can be used to securely cover and seal samples.

PP sealing film,

, which covers a temperature range of -40 to +120°C, is particularly suitable for storage for PCR and ELISA.

Aluminum sealing films are suitable for cold storage due to their good temperature resistance from -80 to +120°C.

## Additional Information

---

No.	BN-2269
Manufacturer (Brand)	Brand
old neolab article no.	7-5237
EAN	4058072214807
Transport temperature	Room temperature
Color	silver
Material	Aluminum
Temp MAX	120 °C
Temp MIN	-80 °C
sterile	No
Format/ Packing	individually packed

