



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® storage adhesive film

Special Price

€34.60 was

€43.60

**plus VAT &
Shipping**

Product Images



Description

The neoLab- sealing adhesive films securely seal your multiwell plates, thus avoiding contamination.

They are suitable for various applications: PCR, qPCR, storage, cell culture, etc..

The adhesive foils are easy to apply and remove again without adhesive residues.

All adhesive foils are free of RNase, DNase and human DNA.

PU: 100 pcs./pack

Additional Information

No.	7-5213
Manufacturer (Brand)	neoLab
EAN	4058072041458
Transport temperature	Room temperature
Color	transparent (milky)
Material	Polyester (PEs)
Temp MAX	110 °C
Temp MIN	-20 °C
sterile	No
Wide	80 mm
Length	80 mm
Dimensions	80 x 130 mm

