



neoLab® optical qPCR adhesive film

Product Images



Description

The neoLab- sealing adhesive films securely seal your multiwell plates, thus avoiding contamination. They are suitable for various applications: PCR, qPCR, storage, cell culture, etc.. The adhesive foils are easy to apply and remove again without adhesive residues. All adhesive foils are free of RNase, DNase and human DNA.

PU: 100 pcs./pack

Additional Information

No.	7-5216
Manufacturer (Brand)	neoLab
EAN	4058072041472
Color	transparent (kristallklar)
Temp MAX	110 °C
Temp MIN	-40 °C
sterile	No
Wide	77 mm
Length	77 mm
Dimensions	77 x 140 mm

