



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® Autosampler Vials 4 ml, clear glass,  
45 x 15 mm, ND13, 100 pcs./pack**

**€27.10  
plus VAT &  
Shipping**

**Product Images**



## Description

---

### **Autosampler vials 4 ml with thread 13-425.**

Made of borosilicate glass type I.

Height 45 mm, diameter 15 mm.

## Additional Information

---

No.	2-7210
Manufacturer (Brand)	neochrom
EAN	4058072009045
Transport temperature	Room temperature
Volume	4 ml
Color bottle	transparent (clear)
Height	45 mm
Floor outside	flat
Floor inside	flat
Thread	13-425 (US)
Neck size	far
DM bottle	15 mm
DM Neck	13 mm
DM Closure	13 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
ND (nominal diameter)	ND13
for the autosampler model	KA

