



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® short thread vials ND9, clear glass,  
with integrated 0.2 ml micro insert, 32 x 11.6 mm,  
with labelling field and filling mark, 100  
pcs./pack**

**Price on  
request  
plus VAT &  
Shipping**



## Description

---

The bottles can be used on all common autosamplers due to their technical geometry.  
The wide opening allows easy filling with viscous substances.

Dimensions: 32 x 11.6 mm

## Additional Information

---

No.	7-0679
Manufacturer (Brand)	neochrom
EAN	4058072029425
Transport temperature	Room temperature
Volume	0.2 ml
Color bottle	transparent (clear)
Height	32 mm
Floor outside	flat
Floor inside	flat
Thread	Short thread
Neck size	far
Locking micro insert	at the bottleneck
DM bottle	11.6 mm
DM Closure	9 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
ND (nominal diameter)	ND9
for the autosampler model	Agilent 1050 Agilent 1090 Agilent 1100 Agilent 1200 Agilent 5880 Agilent 5890 Agilent 6850. 27 Pos. Tray Agilent 6890 Agilent 7673A Agilent 7683A Agilent 79855(A) Agilent CTC GC PAL Agilent CTC HTS/ HTC PAL AI 42 vial tray AI 60 vial tray AI CTC A200S AIM

