



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® thread vial ND10 1.5 ml, clear glass,  
wide opening, 32 x 11.6 mm, 100 pcs./pack**

**€21.70  
plus VAT &  
Shipping**



## Description

---

Standard wide-mouth clear glass sample vial.  
1.5 ml capacity, 10-425 thread.  
Made of 1st hydrolytic grade glass.

## Additional Information

---

No.	7-0725
Manufacturer (Brand)	neochrom
EAN	4058072155209
Transport temperature	Room temperature
Volume	1.5 ml
Color bottle	transparent (clear)
Height	32 mm
Floor outside	flat
Floor inside	flat
DM bottle	11.6 mm
DM Closure	10 mm
Thread	10-425 (US)
Neck size	far
Material bottle	Borosilicate glass (1st hydrolytic class)
ND (nominal diameter)	ND10
for the autosampler model	Alcott 719 AL Alcott 719 D/ D-PCS Alcott 738 Antec-Leyden AS 100 Antek 736 Unisampler Antek 738 Beckmann 501 Beckmann 502/ 502e Beckmann 507/ 507e Cambridge Scientific Instruments 205 Series Cambridge Scientific Instruments 300 Series Cecil Instruments Au

