



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® thread bottles amber glass 1,5ml
32x11,6mm, w.title block + fill mark. Pack a 100
pieces**

**€15.50
plus VAT &
Shipping**



Description

Standard wide-mouth amber glass sample vial.
1.5 ml capacity, thread 10-425.
Made of 1st hydrolytic grade glass.

Additional Information

No.	7-0727
Manufacturer (Brand)	neochrom
EAN	4058072029814
Transport temperature	Room temperature
Volume	1.5 ml
Color bottle	brown
Height	32 mm
Floor outside	flat
Floor inside	flat
Thread	10-425 (US)
Neck size	far
DM bottle	11.6 mm
DM Closure	10 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Alcott 719 AL Alcott 719 D/ D-PCS Alcott 738 Antec-Leyden AS 100 Antek 736 Unisampler Antek 738 Beckmann 501 Beckmann 502/ 502e Beckmann 507/ 507e Cambridge Scientific Instruments 205 Series Cambridge Scientific Instruments 300 Series Cecil Instruments Au

