



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Wheaton sample vials amber glass, 1.5 ml, weight 8 mm, 200 pcs/pack

Special Price

**€53.00** was

**€126.38**

**plus VAT &  
Shipping**

### Product Images



## Description

---

### Sample vials amber glass

High-quality sample vials made of borosilicate glass resistant to thermal shock. Autoclavable.

Packed in trays of 100 pieces each, dust-free.

Amber glass protects light-sensitive samples and maintains purity of contents.

Supplied without screw caps.

## Additional Information

---

No.	WD-0673
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0158
EAN	4058072426873
Transport temperature	Room temperature
Volume	1.5 ml
Color bottle	brown
Height	32 mm
Floor outside	flat
Floor inside	flat
Thread	08-425 (US)
Neck size	narrow
DM bottle	12 mm
DM Neck	8 mm
DM Closure	8 mm
Material bottle	Borosilicate glass (1st hydrolytic class)

