



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Threaded vial w. lid, soda lime glass, 50 x 27 mm, 20ml, DIN 26

€54.90
plus VAT &
Shipping

Product Images



Description

The Ecolab threaded vials with lids are made of soda lime glass and can be used in a variety of ways.

Additional Information

No.	E-0050
Manufacturer (Brand)	neoLab
EAN	4058072082802
Transport temperature	Room temperature
Volume	20 ml
Color bottle	transparent (clear)
Height	50 mm
Floor outside	flat
Floor inside	flat
Thread	DIN 25
Closure / lid properties	closed
DM bottle	27 mm
Color closure/ lid	black
Material bottle	Soda-lime glass (normal glass) 2

