



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Threaded vial w. lid, soda lime glass, 90 x 30 mm, 50 ml, DIN 26

**€70.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The Ecolab threaded vials with lids are made of soda lime glass and can be used in a variety of ways.

## Additional Information

---

No.	E-0053
Manufacturer (Brand)	neoLab
EAN	4058072082833
Transport temperature	Room temperature
Volume	50 ml
Color bottle	transparent (clear)
Height	90 mm
Floor outside	flat
Floor inside	flat
Thread	DIN 25
Closure / lid properties	closed
DM bottle	30 mm
Color closure/ lid	black
Material bottle	Soda-lime glass (normal glass) 2

