



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Vials 2 ml clear glass, cap closed,
septum butyl/PTFE, 100 pcs./pack**

**€109.20
plus VAT &
Shipping**

Product Images



Description

Sample vial made of clear borosilicate glass type I (5.1)
with screw cap, complete with butyl rubber liner with PTFE coating,
temperature resistant up to 110° C.
Volume: 2 ml

Additional Information

No.	2-7230
Manufacturer (Brand)	neoLab
EAN	4058072009083
Transport temperature	Room temperature
Volume	2 ml
Color bottle	transparent (clear)
Height	36 mm
Floor outside	flat
Floor inside	flat
Quantity components	a
Closure / lid properties	closed
DM bottle	12 mm
Color closure/ lid	black
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Polypropylene (PP)
TBST MAX Septum	110 °C

