



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® Vials 30 ml clear glass, cap closed,  
septum butyl/PTFE, 100 pcs./pack**

**€170.40  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Sample vial made of clear borosilicate glass type I (5.1)  
with screw cap, complete with butyl rubber liner with PTFE coating,  
temperature resistant up to 110° C.  
Volume: 30 ml

## Additional Information

---

No.	2-7236
Manufacturer (Brand)	neoLab
EAN	4058072009144
Transport temperature	Room temperature
Volume	30 ml
Color bottle	transparent (clear)
Height	75 mm
Quantity components	a
Floor outside	flat
Floor inside	flat
DM bottle	27 mm
Color closure/ lid	black
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Polypropylene (PP)
TBST MAX Septum	110 °C
Closure / lid properties	closed

