



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton sample vial clear glass, without cap, 24 ml PU 200 pcs. No.9-0157

Special Price

€254.66 was

€314.76

**plus VAT &
Shipping**



Description

High-quality sample vials made of borosilicate glass resistant to temperature changes

Autoclavable

Packed in trays of 100 pieces each, dust-free.

Amber glass protects light-sensitive samples and maintains purity of contents.

Supplied without screw caps.

Additional Information

No.	WD-0707
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0156
EAN	4058072427030
Transport temperature	Room temperature
Volume	24 ml
Color bottle	transparent (clear)
Height	85 mm
Floor outside	flat
Floor inside	flat
Thread	20-400 (US)
Neck size	far
DM bottle	23 mm
DM Closure	20 mm
Material bottle	Borosilicate glass (Wheaton 33. ASTM standard type 1 class A. USP standard type 1)
TBST MAX bottle	505 C°

