



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Wheaton sample vial, amber glass, without cap,  
4 ml PU 200 pcs. No.9-0161**

**€157.66  
plus VAT &  
Shipping**

## Product Images



## Description

---

High-quality sample vials made of borosilicate glass resistant to temperature changes

### **Autoclavable**

Packed in trays of 100 pieces each, dust-free.

Amber glass protects light-sensitive samples and maintains purity of contents.

Supplied without screw caps.

## Additional Information

---

No.	WD-0709
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0160
EAN	4058072427054
Transport temperature	Room temperature
Volume	4 ml
Color bottle	brown
Height	45 mm
Floor outside	flat
Floor inside	flat
Thread	13-425 (US)
Neck size	far
DM bottle	15 mm
DM Closure	13 mm
Material bottle	Borosilicate glass (Wheaton 320. USP standard type 1)
TBST MAX bottle	510 C°

